

**Mycoplasma genitalium In Vitro Transcript, RNA****Analyte Specific Reagent (ASR).**

Analytical and performance characteristics are not established.

**Intended Use**

Analyte Specific Reagent. Analytical and performance characteristics are not established. This analyte specific reagent is intended to be used as a component of a laboratory developed test (LDT) used exclusively by laboratories. Laboratories can develop their own LDTs using one or more Analyte Specific Reagents and/or general laboratory reagents for *in vitro* diagnostic purposes.

**Materials Provided**

Mycoplasma genitalium In Vitro Transcript, RNA (Cat. No. PRD-03410)

Component	Quantity	Concentration	Description
Tube	1 x 2 mL	50,000 copies IVT/mL	<i>In vitro transcript of M. genitalium 16S rRNA.</i>

**Warnings and Precautions**

- A. For Laboratory Use only. Analytical and performance characteristics are not established.
- B. Use routine laboratory precautions. Wear disposable, powderless gloves, protective eye wear, and laboratory coats when handling specimens and kit reagents. Wash hands thoroughly after handling reagents.
- C. Dispose of all materials that have come in contact with reagents according to local, state, and federal regulations.
- D. For information on any hazard and precautionary statements that may be associated with reagents refer to the Safety Data Sheet Library at [www.hologicsds.com](http://www.hologicsds.com).

**Product Specifications****Quality Control**

This analyte specific reagent was manufactured and released in accordance with the Hologic quality control and quality assurance procedures. This reagent has been functionally tested using a system capable of detecting luminescence along with commercially available transcription mediated amplification (TMA) general purpose reagents (GPR). However, it is the responsibility of each laboratory to develop and validate their own tests for clinical diagnostic purposes.

**Storage and Stability**

Please refer to product labels for information on storage and stability.



Hologic, Inc.  
10210 Genetic Center Drive  
San Diego, CA 92121 USA

Customer Support: +1 844 Hologic (+1 844 465 6442)  
customersupport@hologic.com

Technical Support: +1 888 484 4747  
molecularsupport@hologic.com

For more contact information visit [www.hologic.com](http://www.hologic.com).

Hologic and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

All other trademarks that may appear in this package insert are the property of their respective owners.

This product may be covered by one or more U.S. patents identified at [www.hologic.com/patents](http://www.hologic.com/patents).

© 2017 Hologic, Inc. All rights reserved.

AW-16923 Rev. 001  
2017-07